



Reference: PTXBC014940

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	2410002I16Rik; 5830420C16Rik; CFP1; CGBP; HsT2645; PCCX1; PHF18; SPP1; ZCGPC1; hCGBP; CXXC-type zinc finger protein 1; CXXC finger 1 (PHD domain); DNA-binding protein with PHD finger and CXXC domain; PHD finger and CXXC domain-containing protein 1; cpG-bi		

For research use only.

GENE DESCRIPTION	CXXC finger 1 (PHD domain)
GENE ID	30827
ORIGIN SPECIES	Human
ACCESSIONS	BC014940
VECTOR	pDONR223
NCBI SYMBOL	CXXC1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]