



Reference: PTXBC014963

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	pelle; interleukin-1 receptor-associated kinase 1; IRAK-1; Pelle homolog; interleukin 1 receptor associated kinase 1		

For research use only.

GENE DESCRIPTION	interleukin-1 receptor-associated kinase 1
GENE ID	3654
ORIGIN SPECIES	Human
ACCESSIONS	BC014963
VECTOR	pDONR223
NCBI SYMBOL	IRAK1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]