



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

XPA-xeroderma pigmentosum, complementation group A Gene

Reference: PTXBC014965

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	XPA, DNA damage recognition and repair factor; XP1; XPAC; DNA repair protein complementing XP-A cells; excision repair-controlling; mutant xeroderma pigmentosum complementation group A; xeroderma pigmentosum group A-complementing protein; xeroderma pigmen		

For research use only.

GENE DESCRIPTION	xeroderma pigmentosum, complementation group A
GENE ID	7507
ORIGIN SPECIES	Human
ACCESSIONS	BC014965
VECTOR	pDONR223
NCBI SYMBOL	XPA
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr