

LENG1-leukocyte receptor cluster (LRC) member 1 Gene

Reference: PTXBC014986

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	leukocyte receptor cluster member 1; leukocyte receptor cluster (LRC) member 1		
For research use only.			

GENE DESCRIPTION	leukocyte receptor cluster (LRC) member 1
GENE ID	79165
ORIGIN SPECIES	Human
ACCESSIONS	BC014986
VECTOR	pDONR223
NCBI SYMBOL	LENG1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaatatcttgcacgaagagctggcagctccggaacaaggacaatgtcgcccgctgcccgtgacgagcccgagccgggaggaggagaaggagcgtgagcggagggtgctgctgctcagcaaggccgtacagaattcctacggaagaagccagacatcagaa
ctcactgctgagcttgaagcagcagaggcgggagcccgaggttctgcccgtgagacctgttgcgggagctgctgagggaagggaaggagtgatcagaggcaataaagagtagcagggaagaaagcgacaggagaagagaggaagagaagctctggcctctgacatac
ctggccagagtgacgaggacagactcaacccctgtgtacagctaccccgagggcgagggggcccccgcggcccgccagccagatgagaagatcaagagccgtctggaccctctgcccagatgcagaagcatctggggaagaagagacagcacggcgtgatga
aggcagtcgagcagaaggaagagggtctgagaagcagcgacccaaggagcctccatccctggaccagcttcagctgaacgtctgaggagggaagcagctgagaggtctcggcagaggccctgctggccgggtccaaggccggcactacaggagggtcagccg
gaagaagacgagacggatgaccggcgggcggtacaaactccaattcaaccccgagctggcccgccggcccgccagcaggacctacactactactga