

RAB23-RAB23, member RAS oncogene family Gene

Reference: PTXBC015021

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RAB23, member RAS oncogene family; HSPC137; ras-related protein Rab-23; RAB family small GTP binding protein RAB 23		

For research use only.

GENE DESCRIPTION	RAB23, member RAS oncogene family
GENE ID	51715
ORIGIN SPECIES	Human
ACCESSIONS	BC015021
VECTOR	pDONR223
NCBI SYMBOL	RAB23
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgttgaggagaagatggaagtcgcatagagatgggtgttgaggatggagcagttggaaaatcaagtagattcagcgatattgcaaggcattttacaaaagactacaagaaaccattggagtgatttttggagcgacaattcaagttaatgatgaagatgcagactaatgttat
gggacactgcaggtcaggaggaatttgatgcaattacaaaggcctactatcgaggagcccaggtctgtgtgctcgtgttctctaccacagatagggaatcttttgaagcagtttccagttggagagagaaaagtagccgaagtgggagatataccaactgtactgtgcaaaacaagattg
atcttctggatgattctgtataaagaatgagggaagctgaggcactggcaaaagggttaaggttaagattctacagaacatcagtgaaagaagatcctaaatgtgaatgaagtgttttaagtatttggctgaaaaataccttcagaaactcaaacacaataagctgaggatccagaactaacgcatt
caagtagtaacaagattggtgtctttaatacatctggtggaagtcactccggtcagaattcaggtaccctcaatggtggagatgtcatcaatcttagaccaacaacaaggaaccaagaaaaacagaaatccttttagcagctgtagcataccctaa