

**UBR7-ubiquitin protein ligase E3 component n-recognin 7 (putative) Gene**

Reference: PTXBC015046

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C14orf130; N-recognin-7; ubiquitin protein ligase E3 component n-recognin 7 (putative)		

**For research use only.**

---

<b>GENE DESCRIPTION</b>	ubiquitin protein ligase E3 component n-recognin 7 (putative)
<b>GENE ID</b>	55148
<b>ORIGIN SPECIES</b>	Human
<b>ACCESSIONS</b>	BC015046
<b>VECTOR</b>	pDONR223
<b>NCBI SYMBOL</b>	UBR7
<b>LENGTH</b>	ProteoGenix bp
<b>SEQUENCE PRIMERS</b>	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

**ORF (Only for reference, product provided without STOP codon)**

atgatccagtcgtagtctgtgaagactgttcctcatggaagcatctgtggtccattcccctgagagtgaggatcttcaggagatgtagccaggcctcatgaaactgttctttttgtggccttagctgcacaattggcagtaaccccaaaatccactgaggatgagattgtgcgg  
aacattgatggaatagtgatcaggaaagttaacaaactgaaaatggagagcatcaagatagtagtaccctcaagaggatgtccagaacaggaaaggatgatgtccgggaggttaaagttagagcagaacagtgaacctgtgccgctctagttctgaatctgatctccagacagtgttaa  
gaatgaaagcctcaacgcagaatcaaaatctggctgcaaaactcaggagcttaagctaagcagcttataaagaaagacactgccacattggcccctgaaactggcgtagcaagttgtgtacctgccaagactgtatgaaatgtatggagatctagatcttattctgacagatgaatac  
gacacagttctggcttatgaaacaaagggaagatgccagccactgacagagcgcattccctaatggatacccttagcagcatgaatagagtcagcaagtggaactcattgtgaatacaatgattgaaactgaactaaagactatcagaagattgctgatgaaagcagcgtt  
gttaagagagagacattcagcagttcttgaagatttcagtcaaaaaagagaagagagtgatggatgcagtattactgcagctag