

**MOBKL2A-MOB1, Mps One Binder kinase activator-like 2A (yeast) Gene**

Reference: PTXBC015049

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MOBKL2A; MOB-LAK; MOB1C; moblak; MOB kinase activator 3A; MOB LAK; MOB1, Mps One Binder kinase activator-like 2A; mob1 homolog 2A; mps one binder kinase activator-like 2A		

**For research use only.**

---

<b>GENE DESCRIPTION</b>	MOB1, Mps One Binder kinase activator-like 2A (yeast)
<b>GENE ID</b>	126308
<b>ORIGIN SPECIES</b>	Human
<b>ACCESSIONS</b>	BC015049
<b>VECTOR</b>	pDONR223
<b>NCBI SYMBOL</b>	MOBKL2A
<b>LENGTH</b>	ProteoGenix bp
<b>SEQUENCE PRIMERS</b>	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

**ORF (Only for reference, product provided without STOP codon)**

atgtccaacccttctgaagcaagtctcaacaaggacaagacattccgcccaagcgaagtttgagccggcaccagcgttcgagctgcacaagaaggcgcaggcgtcgtgaacgccgggctggacctgcggctggccgtgcagctgccccggcgagacctgaacga  
ctgggtggctgttcacgtgggtggaactctttaaccgcgtcaacctcatctacggaccatcagcagcggctgcacggagcagctcgcgccgtatgtcggggggcccaagtatgagtagcgcctggcaggatgagcataaagtccggaagcccacggcactgtccgcgccaggtacat  
ggacctgctgatggactggatcggggcgagatcaacaacgaggacctctccccaccaacgttggcactccgttccaagaactcctgcagacgggtgcggaagatcctgcggcgtgtccgcgtgttcgtgcacgtctacatccaccattgaccgcacgcgagatggctcc  
gaggcccacgtgaacacctgctacaagcactctactatttctcaaggagttcggcctcatcgacaccaaggagctggagccactgaaagaaatgaccgcccgatgtgccactga