



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

OIP5-Opa interacting protein 5 Gene

Reference: PTXBC015050

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

5730547N13Rik; CT86; LINT-25; MIS18B; MIS18beta; hMIS18beta; protein Mis18-beta; LAP2alpha interactor-25; MIS18 kinetochore protein homolog B; cancer/testis antigen 86; opa-interacting protein 5; Opa interacting protein 5

For research use only.

GENE DESCRIPTION	Opa interacting protein 5
GENE ID	11339
ORIGIN SPECIES	Human
ACCESSIONS	BC015050
VECTOR	pDONR223
NCBI SYMBOL	OIP5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr