

CBFA2T2-core-binding factor, runt domain, alpha subunit 2; translocated to, 2 Gene

Reference: PTXBC015066

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein CBFA2T2; EHT; MTGR1; ZMYND3; p85; ETO homolog on chromosome 20; ETO homologous on chromosome 20; MTG8-like protein; MTG8-related protein 1; core-binding factor, runt domain, alpha subunit 2; translocated to, 2; myeloid translocation gene-related p		

For research use only.

GENE DESCRIPTION	core-binding factor, runt domain, alpha subunit 2; translocated to, 2
GENE ID	9139
ORIGIN SPECIES	Human
ACCESSIONS	BC015066
VECTOR	pDONR223
NCBI SYMBOL	CBFA2T2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

```
atgcatggatcgctgtggaagtgaagatacagtcagatctcactcccacatgccacccctcccaccaataaacctggaggaccgagccagtgctctcactctactgcatlaagcaatggccatcaaccattctcctctaccctgaatgggcccacaccgccacagagattc
agcaatggtcctgctcctccacatcctcactcacaatacgaatggccagccactgtggtgctcgacaactcagcaagtgaaacgcttctaccactcgaacagttggcaatgacatccctgagatggggagaaggtgcgggactctgttcttgcactgtgaactcaaca
gtgacaattgaggaaatccactgtaagctccaagaagccacaacttcccctcgtcctttgtgattccatttcaaggccaacctgccctcgtcagcgggaaactgctgactgcctcggcgccaagcagaccatcccagtaacctggctcagcacgaacacctctgctcaac
acaagcattgcatgcctgctgactcagagttgctcaggtgcacggaatgggaagggccagtcagagagctcataatacactaaaatctatgggagtagcttcagacatgagtattttgggaactctcagatgcctgcacaccggtaagaactctacttagag
aatttatggctgactgtccagaagcctcaggagtagtctgcttttcctga
```