

## NT5DC1-5'-nucleotidase domain containing 1 Gene

Reference: PTXBC015138

|                |                                                                                                                                                           |                     |                |
|----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                       | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                      | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | C6orf200; LP2642; NT5C2L1; 5'-nucleotidase domain-containing protein 1; 5'-nucleotidase, cytosolic II-like 1 protein; 5'-nucleotidase domain containing 1 |                     |                |

**For research use only.**

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | 5'-nucleotidase domain containing 1                                                               |
| <b>GENE ID</b>          | 221294                                                                                            |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC015138                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | NT5DC1                                                                                            |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon)

atgtcagcttttaaggaaaactgtggaatatatttccagaataaaaaagagatcagcagatattacatagttgtcctgaatctgtgaaaaatggcttcgacagctaaagaatgctgggaaaattcttctgtaattaccagttctcacagtgttagcttctgcgaatatattctgg  
gaatgatittacagacctttttgacattgtgattacaaatgcattgaagcctgtttcttctccacttaccaagtcagagacctttccggacactcgagaatgatgaggagcaggaggcactgccatctctggataaacctggctgtactcccaagggaacgctgtccacctctatgaactctg  
aagaaaatgactggcaaacctgaaccaaggtgtttattttgtgacagcatgcatcagatatttccagctcgtcactatagtaattgggagacagctcctcatctgggaagaactcagaggggatgaaggcagaggagtcagaggcctcagagagaaag  
gaaaatagagggacaaaagcaaacctttaaatactcatctaaaaatggggctctttttattgattcagttttgggactggaaaatacagaagactcctgtttatcatggtcttgaagaatcagctactacagcactattgcaattccaagtattgaagcaatcgagaattacctctg  
gactacaaattacaagattctctcaagcaattcaaaaacagctggctactatccaaatcctccactggcttatcaagtgatgagacactgataccaataa