

C13orf27-chromosome 13 open reading frame 27 Gene

Reference: PTXBC015148

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C13orf27; testis-expressed sequence 30 protein; testis expressed 30		
For research use only.			

GENE DESCRIPTION	chromosome 13 open reading frame 27
GENE ID	93081
ORIGIN SPECIES	Human
ACCESSIONS	BC015148
VECTOR	pDONR223
NCBI SYMBOL	C13orf27
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaattctctcattgatgtcactggcatcccatcttgcatctcatgggttttctgctgagatttacctgtaaaaggccttaattgtacatagaattaaggcgtataaatcagttttgaattacctgaagacatcaggagaatacaaaactgcaggtgttttctggaggtcgttcaatgggtcaa
gagcagctgcttctgaatgtgtcacattgagccagatgatgtgatgtttgtcggtgtcattgtatttctaccactgcaccatccaaagcagcagcataaactcagagatgaagacctcttctgtttaaagagcctgtactgtttgtcaggtcagcagatgaaatgtgtgaaaaga
acttgttgagaaagtggcacagaaaatgcaagctcccataaaatccactgattgagaaggcaaatcattccatggcagtgaaaggacggtcgacaaatgatgtttcaagaataatacacagattttgtttgatccaagaataactgaaatggacaagaatgtcattag