

C3orf49-chromosome 3 open reading frame 49 Gene

Reference: PTXBC015210

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	chromosome 3 open reading frame 49		

For research use only.

GENE DESCRIPTION	chromosome 3 open reading frame 49
GENE ID	132200
ORIGIN SPECIES	Human
ACCESSIONS	BC015210
VECTOR	pDONR223
NCBI SYMBOL	C3orf49
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcccaacctcagctgtatctaccagaaccctttaagattgctacagaaaagttggacagtgccgaagattccagcaactcaagaagaaaatggctcattcaaaaggaagggatagaaagatggcatagagctgtgtctaccaacttactaaaaaaatgtattgtccctaaagag
gaatcatcagtgatagtgacattgggatttcataagccagcaaaatcagaaaagtaattgaagacgaaagtgaagactgctttgggaggtatgctgtctacaagtatagaagcaagccagcatgtgccagccaaggggctccactgaccataaagaagccctcctgtcaaacacgc
agagccttcctcctaggtgtcaaggagtttcatctcccaagttatttacagcaaaaatgagaaagctctcagagaatgcgactatacaactgtatgtgtagagccagagacagaggagataaccagggaacacactcctcgggagagaaacccaagcggttatctgtgacat
ctcttcctcagagctgcaaaaggggcatattcaccaaaaagagaccacatttccagcacttaaaaaaagaagcgtggcatggaaaacatccttgaaaatcagatttgacagtaggaaagcttcaatgcaggtggatgacctcatagaacagtgactgataaatccatgaagctac
tggtcaaaagacatgctgagctcaacagtgtagtttctggggatgaaattctcagctcttctaaacagttccagaggatattccaagagaacctgaggagataaaataaaatgaccacaaaaggacctggagacagctga