

KIF9-kinesin family member 9 Gene

Reference: PTXBC015311

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	kinesin-like protein KIF9; kinesin protein 9; kinesin family member 9		
For research use only.			

GENE DESCRIPTION	kinesin family member 9
GENE ID	64147
ORIGIN SPECIES	Human
ACCESSIONS	BC015311
VECTOR	pDONR223
NCBI SYMBOL	KIF9
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggtactaggaagagttcatgtcattgtccgtgtcaaacccaccgatgacttgcctcatgaaatgatcagatcggagatgacaaaagaagcattgatacttaaaaaaacattcggagaggagttgtcaataaccaacagacagactggcgttaagttggatggagttcttcac
gatgcctcccaggacttggtttatgagacagttgcaaggatgtgtttctcaggccctcgaaggctataatggttaacttaattaataa