

C7orf31-chromosome 7 open reading frame 31 Gene

Reference: PTXBC015397

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C7orf31; chromosome 7 open reading frame 31		
For research use only.			

GENE DESCRIPTION	chromosome 7 open reading frame 31
GENE ID	136895
ORIGIN SPECIES	Human
ACCESSIONS	BC015397
VECTOR	pDONR223
NCBI SYMBOL	C7orf31
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaccccttgaaactggatgattaccatgaaaagatggtagcagagtaacacgaaagataggattgacccagagccgcaagaaaaattccacctgtcttcaaacctccagaccgttgaaggacgaattgccgattaattcagaaccgacgttctctggaggctatagtcagca
aagaccacgttctgtccagattgtactcctagagttgtgtaattttcatacctttgtacccagttctaaagaaatggtagcacttagtgataacataccagccgtgtgacccataaaaaccaggatattgaagagaagatcatagaagaacagagcttctgtctacatgaactccatctt
gttaccacaacaaagacctaactagcatttatgacataaagccatttccaaaaatcacagatactaaaagacagaagatttactggagacagcagtcactaaaaacccaacccacaccttactgtaaacagacacactggattcctatgaaatcttaaatccccctacgtatcagtata
atatgtgtccagaccctgttagccttagtaaacctagtgttttacaataaacaagacacgggaagcctttcactttagaacatttttaagtaagccagaagaagagttgttcttgaaatggaaacaatgaagaaacaagacgttcttctgttggtgattcctagagctggagtaccaaactcag
accaacctgtctggagcttaagaactcttttcaaaaactggtgcacaaaagcgtttccataaatcaattctagaagaccataaagacctcaggataatgagcatcgggatgaagcaccaattctatggccataattctattttataattga