

ROPN1B-ropporin, rophilin associated protein 1B Gene

Reference: PTXBC015413

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ropporin-1B; AKAP-binding sperm protein ropporin; rophilin-associated protein 1B; ropporin, rophilin associated protein 1B; testis tissue sperm-binding protein Li 83P; rophilin associated tail protein 1B		

For research use only.

GENE DESCRIPTION	ropporin, rophilin associated protein 1B
GENE ID	152015
ORIGIN SPECIES	Human
ACCESSIONS	BC015413
VECTOR	pDONR223
NCBI SYMBOL	ROPN1B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtggaagtggatctccacagatctgttaatagtgatgaatgtgggtcgttcacggaggagatcagtggtctgaagtttttagcccttgctgcagcgtctgggagttactattacaaaactctcaagatgtgtgaggtcttcatgtgaccacaatgggtgggtgcccc
gaatcccatcagcaccttcagtttctacacgtatattgccgaagtggatggggagatctgtgcatcacatgtcagcagatgctaaactacattgaacaggaagtaattggtctgatggtttaatcacgggtgaatgactttacccaaacccagggtttgctggagtaa