

CALB2-calbindin 2 Gene

Reference: PTXBC015484

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CAB29; CAL2; 29 kDa calbindin; calbindin 2, (29kD, calretinin); calbindin D29K; testicular secretory protein Li 8; calbindin 2		

For research use only.

GENE DESCRIPTION	calbindin 2
GENE ID	794
ORIGIN SPECIES	Human
ACCESSIONS	BC015484
VECTOR	pDONR223
NCBI SYMBOL	CALB2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctggcccgagcagcagccccctacctgcacctggccgagctgacggcgtcccagttcctggaaatatggaagcactttgacgcagacggaaatgggtatattgaaggtaaagagctagaaaacttttcaagagctggagaaggcaaggaaaggctctggcatgatgtcaaa
gagtgacaactttggagaaaaatgaaggagttcatgcagaagtatgataaaaactcagatgggaaaaatcgagatggcagagctggcgcagatcctgccaaaccgaagagaacttccttctgtgcttcagcagcagtggtccagcacgagtttatggagcttggcggaagtacg
acacagacaggagtggtctacatcgaagccaatgagctcaagggtattcctgtcagacctgctgaagaaggcgaaccggccgtacgatgagcccaagctccaggaatacacccaaccatactacggatgtttgacttgaacggggatggcaaatgggcctctcagagatgtcccgactc
ctgcctgtccaggaacttcctgttaaatctcaggccatgaagctgacctcagaggagttaacgcgatcttcacattttacgacaaggatagaagcggctacattgacgacatgagctggatgcccttttgaaggatctgtacgagaaaaacaaaaggaaatgaatattcaacagctca
ccaactacagaaagagcgtcatgtccttgcagaggcagggaagctctaccgcaaggacctggagattgtctctgcagcagcccccatgtaa