

IFI27-interferon, alpha-inducible protein 27 Gene

Reference: PTXBC015492

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FAM14D; ISG12; ISG12A; P27; interferon alpha-inducible protein 27, mitochondrial; 2310061N23Rik; ISG12(a); interferon alpha-induced 11.5 kDa protein; interferon-stimulated gene 12a protein; interferon alpha inducible protein 27		

For research use only.

GENE DESCRIPTION	interferon, alpha-inducible protein 27
GENE ID	3429
ORIGIN SPECIES	Human
ACCESSIONS	BC015492
VECTOR	pDONR223
NCBI SYMBOL	IFI27
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaggcctctgctcactcagcagtgaccagtggtggccaaagtggcagggtggcctctggcctgcccgtagttttgccctggccaggattgctacagttgtgattggaggagttgtggctgtgccatgggtgctcagtgccatgggcttcactgcggcggaatgcctcgtcctccatagcagccaagatgatgtccgcggcgccatggccaatgggggtggagttgctcgggcagccttggctactctgagtcactgggagcaactggactctccgattgaccaagttcatcctgggctccattgggtctgccattgcggctgtcattgcaggttctactag