

MSX2-msh homeobox 2 Gene

Reference: PTXBC015509

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CRS2; FPP; HOX8; MSH; PFM; PFM1; homeobox protein MSX-2; homeobox protein Hox-8; msh homeo box 2; msh homeobox homolog 2; msh homeobox 2		

For research use only.

GENE DESCRIPTION	msh homeobox 2
GENE ID	4488
ORIGIN SPECIES	Human
ACCESSIONS	BC015509
VECTOR	pDONR223
NCBI SYMBOL	MSX2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcttctccgtccaaaggcaatgactgttttgcggcgacgaggagggccagcagtggtggccggaccagggccggggcctggggcgccgagggggccgaggaggcgccgctcaaggtctccagcctgcccttcagcgtggaggcgctcatgtccgacaagaagccg
cccaaggaggcgctccccgtgcccggccgaaagcgctcggccggggccacccctgcggccactgctgctgtcggggcgacggcgctcgggaagcgacagccccggcgctggtgaagcccttcgagaccgctcggtaagtcggaaaattcagaagatggagcggcggtgat
gcaggaaccggcgccgatatcgcggccgcaagacatagagccctaccacctgcacccctgaggaacacacaagaccaatcggaagccgacgcccctttaccacatcccagctcctcgccttggagcgcaaggttcggtcagaacagctacctctccattgcagagcgtgcagagttc
tcagctctctgaacctcacagagaccaggtcaaaatctggttcagaaacgaaggccaaaggcgaaaagactgcaggaggcagaactggaagctgaaaatggctgcaaacctatgctgccctcagctcagctcccttccccatcagctcggccctgcaggcagcgtccata
tatggagcatcctaccggttcacatagactgtgcttccatcccgctgtgggactctatgccacgccagtggtgatggcatgtaccacctgtcctaa