

RGS1-regulator of G-protein signaling 1 Gene

Reference: PTXBC015510

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	1R20; BL34; HEL-S-87; IER1; IR20; regulator of G-protein signaling 1; B-cell activation protein BL34; early response protein 1R20; epididymis secretory protein Li 87; immediate-early response 1, B-cell specific; regulator of G-protein signalling 1		

For research use only.

GENE DESCRIPTION	regulator of G-protein signaling 1
GENE ID	5996
ORIGIN SPECIES	Human
ACCESSIONS	BC015510
VECTOR	pDONR223
NCBI SYMBOL	RGS1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccagggaatgttcttctgctaaccacaaagggaattgaagggaaccactcattcactctagacgacaaaatgcaaaaaggaggccaaagacttttggatggatgaaagcatacctgagatctatgatccacatctggaatctggaatgaaatcttccaagtccaaggatgactttctgctgctgaagtaagcaatgggtctcaatctctggaataactcttggcaaccaactggtcaaatgtcttggaaagtttcctaaagctgaattcagtgaggagaataattgagttctgctggcttgtgaagactataagaaaacagagctctgatctttgcctgtaaagcagaagagatatataagcatttggcattcagatgctgctaacaatcaataatgacttccgcactcgagaatctacagccaagaagattaaagcaccacccccacgtgtttgatgaagcacaaaaagtcataatactcttatggaaaaggactcttatccagggttcctcaatcagatatttacttaaatcttctaaatgacctgcaggctaatagcctaaagtga