

RNASE4-ribonuclease, RNase A family, 4 Gene

Reference: PTXBC015520

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RAB1; RNS4; ribonuclease 4; RNase 4; ribonuclease A B1; ribonuclease, RNase A family, 4; ribonuclease A family member 4		

For research use only.

GENE DESCRIPTION	ribonuclease, RNase A family, 4
GENE ID	6038
ORIGIN SPECIES	Human
ACCESSIONS	BC015520
VECTOR	pDONR223
NCBI SYMBOL	RNASE4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctctgcagaggaccattcattgcttctgttttgcctgacccctgtgggctggggctgggtccagccctctatggccaggatggcatgtaccagcgattctgcggcaacacgtgcaccctgaggagacaggtggcagtgatgctactgcaactgatgacaaagacggaa
gatgactttgtatcactgcaagcgcttcaacacctcatccatgaagatatctggaacattcgtagtatctgcagcaccaccaatccaatgcaagaacggcaagatgaactgcatgaggggtgtagtgaaggtcacagattgcagggaacaggaagttccagggcacccaactgcagat
atcgggcatagcgagcactagacgtgtgtcattgcctgtgagggtaaccacaggtgcctgtgcactttgacggtag