

BEX2-brain expressed X-linked 2 Gene

Reference: PTXBC015522

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein BEX2; BEX1; DJ79P11.1; brain-expressed X-linked protein 2; hBex2; brain expressed X-linked 2		

For research use only.

GENE DESCRIPTION	brain expressed X-linked 2
GENE ID	84707
ORIGIN SPECIES	Human
ACCESSIONS	BC015522
VECTOR	pDONR223
NCBI SYMBOL	BEX2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagtcacaaagaggacgagcgttaacaatctcatcgtgaaatgtcaaccaggaaatgatgaaaaagatgaaaaggagcaagttgctaataaaggggagcccttggccctaccttgaatgttagtgaatactgtgtgcctagaggaaaccgtaggcggttccgcgttaggcagc
ccatcctgcagtatagatgggacataatgcataggcttgagagccacaggcaaggatgagagagagaatggaaggattggggaggaggtgagacagctgatgaaaagctgagggaaaagcagttgagtcatagtttgcgggcagtcagcactgatccctcaccatgacc
atcacgatgagtttgccttatgccctga