

KLK6-kallikrein-related peptidase 6 Gene

Reference: PTXBC015525

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	Bssp; Klk7; PRSS18; PRSS9; SP59; hK6; kallikrein-6; kallikrein 6; neurosin; protease M; serine protease 18; serine protease 9; zyme; kallikrein related peptidase 6		

For research use only.

GENE DESCRIPTION	kallikrein-related peptidase 6
GENE ID	5653
ORIGIN SPECIES	Human
ACCESSIONS	BC015525
VECTOR	pDONR223
NCBI SYMBOL	KLK6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

```
atgaagaagctgatgtgtgctgagctgattgctgcagcctgggcagaggagcagaataaagttggtcatggcgaccctgcgacaagacatcaccctaccagctgccctctacacctgggccacttgctctgtgtgggtcttattccatcactgtgggtcctcacagctgc
ccactgcataaaaccgaatctcaggtctctcctgggaaagcataaccttcggcaaaaggagagttccaggagcagagttctgttccggctgtgatccacctgactatgatgccagccatgaccaggacatcatgctgttgcctggcacgccagccaaactctgaaactc
atccagcccttccctggagaggactgctcagccaaccaccagctgccacatcctggcctgggcaagacagcagatggtgattccctgacaccatccagtgatatacaccactggctcccgtgaggaggtgagcatgcctaccctggccagatcaccagaaactgtgtg
gtgctgggatgagaagtacgggaaggattcctgccaggtgattctggggctccgtggtatgtggagaccacctccgagcctgtgtcatgggtaacatcccctgtggatcaaaaggagaagccaggagtctacaccaactctgcagatacagaactggatccaaaaaccattc
aggccaagtga
```