



Reference: PTXBC015526

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	EMA; GA2; MADD; electron transfer flavoprotein subunit alpha, mitochondrial; alpha-ETF; electron transfer flavoprotein, alpha polypeptide; glutaric aciduria II; electron transfer flavoprotein alpha subunit		

For research use only.

GENE DESCRIPTION	electron-transfer-flavoprotein, alpha polypeptide
GENE ID	2108
ORIGIN SPECIES	Human
ACCESSIONS	BC015526
VECTOR	pDONR223
NCBI SYMBOL	ETFA
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]