

AKAP9-A kinase (PRKA) anchor protein (yotiao) 9 Gene

Reference: PTXBC015533

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AKAP9-BRAF fusion protein; AKAP-9; AKAP350; AKAP450; CG-NAP; HYPERION; LQT11; MU-RMS-40.16A; PPP1R45; PRKA9; YOTIAO; A-kinase anchor protein 9; A kinase (PRKA) anchor protein (yotiao) 9; A kinase (PRKA) anchor protein 9; A-kinase anchor protein 350 kDa; A		

For research use only.

GENE DESCRIPTION	A kinase (PRKA) anchor protein (yotiao) 9
GENE ID	10142
ORIGIN SPECIES	Human
ACCESSIONS	BC015533
VECTOR	pDONR223
NCBI SYMBOL	AKAP9
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaggacgaggagagacagaagaagctggaggccggcaagccaagcttgcagtttcgacaagaaaaagctcagtcggatgggcagagtcctccaagaagcagaaaaaaagagaaaaacgtcaagcagtaaacatgatgtgtcagcacacacatgattgaatattgatcaat
cacagtgtaataaataacataaatagttctcagagagtagaatacaactgtgattcctgaatcacaataatgagaactctacatagtgagagaataaccagtcagcagggtctctctgtggaactggaaagtgaatttcaaccacagcagatgactgcagttcagaggtaaatggttgc
agtttttgatgagaacaggaaagcctacaatttataagggaagaagaatttgggttgatgattctattctgaacaaggagcacaagacagtcgactcatctagagatgatgaaagtgaattggctgggaagcagcatgagattgaagagctaaacagagagctggaagaaatgag
ggttacctatgggactgaaggactgcagcagttacaagaatttgaagctgccattaaacaagaagatggcattataaccagctcactgctaattacaacaagcaagaagagaaaaggatgagacaatgagagaatttttagattgacagaacagagtcacaaattacagattcaatttcag
caattacagctagtgaaactctgagaacagcactcatagtagcacagctgcagacttactacaagccaacaacagatcctcactcatcaacagcagcttgaaagaacaagaccacttattagaagattatcagaaaaagaagaagacttcacaaatgcaaatgatttcttgcagagaaa
attaaagtatatgaatggtatggtttatttaa