

HBE1-hemoglobin, epsilon 1 Gene

Reference: PTXBC015537

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HBE; hemoglobin subunit epsilon; epsilon globin; hemoglobin epsilon chain; hemoglobin, epsilon 1; hemoglobin subunit epsilon 1		

For research use only.

GENE DESCRIPTION	hemoglobin, epsilon 1
GENE ID	3046
ORIGIN SPECIES	Human
ACCESSIONS	BC015537
VECTOR	pDONR223
NCBI SYMBOL	HBE1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtgcatttactgctgaggagaaggctgccgtcactagcctgtggagcaagatgaatgtggaaggctggaggtgaagcctgggcagactcctcgtgtttaccctggaccagagatgtttgacagctttgaaacctgtcgtctcctctgccatcctgggaacccaagggtca
aggcccatggcaagaaggctgacttctttggagatgctataaaacatggacaacctcaagccgcctttgctaagctgagtgagctgcaactgtgacaagctgcatgtggatcctgagaactcaagctcctgggtaacgtgatggtgattattctggtactacatttgcaaggagttc
accctgaagtgcaggctgcctggcagaagctgggtgtctgtcgtccattgccctggcccataagtaaccactga