



Reference: PTXBC015542

|                               |                                                                             |                     |                |
|-------------------------------|-----------------------------------------------------------------------------|---------------------|----------------|
| Product type:                 | DNA                                                                         | Storage Condition:  | -20°C          |
| Size/Quantity:                | 2 ug                                                                        | Delivery Condition: | +4°C, Blue Ice |
| Synonyms                      | HVED; NECL5; Necl-5; PVS; TAGE4; poliovirus receptor; nectin-like protein 5 |                     |                |
| <b>For research use only.</b> |                                                                             |                     |                |

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | poliovirus receptor                                                                               |
| <b>GENE ID</b>          | 5817                                                                                              |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC015542                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | PVR                                                                                               |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

[illegible]