

HIST2H3A-histone cluster 2, H3a Gene

Reference: PTXBC015544

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	H3/n; H3/o; histone H3.2; histone 2, H3a; histone H3/o; histone cluster 2, H3a; histone cluster 2 H3 family member a		

For research use only.

GENE DESCRIPTION	histone cluster 2, H3a
GENE ID	333932
ORIGIN SPECIES	Human
ACCESSIONS	BC015544
VECTOR	pDONR223
NCBI SYMBOL	HIST2H3A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcccgtactaagcagactctcgcaagtcgaccggcggaagggcccgaggaagcagctggccacaaaggcgcccgcaagagcgccggccacggcggggtgaagaagccgcaccgtaccggccggcaccgtagccctgcgggagatccggcgtaccagaa
gtccacggagctgctgatccgcaagctgccctccagcggctgtacgcgagatcgcgaggactttaagacggacctgcgctccagagctcggcgtgatggcgtgcaggaggccagcgagccctacctgtggggctgttcgaagacacgaacctgtgcgcatccacgcaa
gcgctgacattatgcccaaggacatccagctggcccgccatccgtggagagcgggcttaa