



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

AKT1S1-AKT1 substrate 1 (proline-rich) Gene

Reference: PTXBC015562

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	Lobe; PRAS40; proline-rich AKT1 substrate 1; 40 kDa proline-rich AKT substrate; AKT1 substrate 1 (proline rich); AKT1 substrate 1		

For research use only.

GENE DESCRIPTION	AKT1 substrate 1 (proline-rich)
GENE ID	84335
ORIGIN SPECIES	Human
ACCESSIONS	BC015562
VECTOR	pDONR223
NCBI SYMBOL	AKT1S1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]