



Reference: PTXBC015581

|                |                                                                                                                                                                                                                                           |                     |                |
|----------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                                                                                       | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                                                                                      | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | FIGLER6; SALM3; SALM3; leucine-rich repeat and fibronectin type-III domain-containing protein 4; fibronectin type III, immunoglobulin and leucine rich repeat domains 6; leucine rich repeat and fibronectin type III domain containing 4 |                     |                |

**For research use only.**

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | leucine rich repeat and fibronectin type III domain containing 4                                  |
| <b>GENE ID</b>          | 78999                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC015581                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | LRFN4                                                                                             |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

[illegible]