

CPNE3-copine III Gene

Reference: PTXBC015734

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CPN3; PRO1071; copine-3; copine III; copine 3		

For research use only.

GENE DESCRIPTION	copine III
GENE ID	8895
ORIGIN SPECIES	Human
ACCESSIONS	BC015734
VECTOR	pDONR223
NCBI SYMBOL	CPNE3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcttggaaacagcatagatatgtgctgtggttttcagaatttctcttttaacacaagaagccttttaaaaaatgacttacacatatctcatgtacagtaaaacagacagaagtgtgcttatctgtttga