

C1orf124-chromosome 1 open reading frame 124 Gene

Reference: PTXBC015740

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C1orf124; DVC1; PRO4323; spartan; sprT-like domain-containing protein Spartan; DNA damage protein targeting VCP; DNA damage-targeting VCP (p97) adaptor; zinc finger RAD18 domain-containing protein; SprT-like N-terminal domain		

For research use only.

GENE DESCRIPTION	chromosome 1 open reading frame 124
GENE ID	83932
ORIGIN SPECIES	Human
ACCESSIONS	BC015740
VECTOR	pDONR223
NCBI SYMBOL	C1orf124
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggatgatgacttgatgtggcactcggcctcagaggagtggaactgcaggaggcggagcgcgatcatgccaggagtcctgtcgtcgtgacgcgtcgtggagttggtggacccacacgggactgcaggcactgtttgtcagtttaacaccaattctctggggccagct
ggagccgtcgaggtgaagtggagcgtgcgaatgaccctgtgtcctggatagcagctatgaagggaagggtggaatgtgtccatccgtctcagcgaacccctttgaggtgagcccaagaaagatctgttagagccctcctcgtcatgaatgatacatgcctatttattgtcactaat
aacgacaagaccgagaaaggcgtgtccagaattttgtaacatatgcacgcacacagcctgactggagccaataaacggtataaccatacttttcacgatgaggtggtgagtagtgcggcgacactggtggcgtcgaatgggcccgtgccagcacaggccaccgtattacggctatg
tcaaacgagctactaacagggaacccctcgtcctcatgactattgtggcgtgagcaccagaaaacctgtggagcacttacataaaatcaagggaaccagagaattactcaaaaaaggcaaggaaaggcaaaactaggaaaggaaaccagattggccgcagagaataaaggtaaccttc
gtgtaa