

LAGE3-L antigen family, member 3 Gene

Reference: PTXBC015744

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	EKC/KEOPS complex subunit LAGE3; CVG5; DXS9879E; DXS9951E; ESO3; ITBA2; Pcc1; protein ESO-3; L antigen family member 3		

For research use only.

GENE DESCRIPTION	L antigen family, member 3
GENE ID	8270
ORIGIN SPECIES	Human
ACCESSIONS	BC015744
VECTOR	pDONR223
NCBI SYMBOL	LAGE3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctgggacgctgatcagacgagcggaggcgtgacggcgggatggccgggtggccacagctccgctggcggcgtggacacagccgagctccggcgggtggagctccccagcgcacgcgcccaggtccgggcagagacgcccgtctgctggccaggggtcac
gaatgctggccgcacatattcaccctcagcgtgcttccgacccccctggaggcggaaatgccatgggtccctggcaccagatccgagccccacaaagggtgttgggaaggatcacaagtgaftggcagatcctgctgctccctggaaagctgaagactgctgctgctc
cgaattccgcatcaacttcttgaccagcttccctggtgctgacatgacagccttggcggcccccgttccccgctaa