

NKG7-natural killer cell group 7 sequence Gene

Reference: PTXBC015759

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein NKG7; GIG1; GMP-17; p15-TIA-1; G-CSF-induced gene 1 protein; GIG-1 protein; granule membrane protein 17; granule membrane protein of 17 kDa; natural killer cell group 7 sequence; natural killer cell protein 7; natural killer cell granule protein 7		

For research use only.

GENE DESCRIPTION	natural killer cell group 7 sequence
GENE ID	4818
ORIGIN SPECIES	Human
ACCESSIONS	BC015759
VECTOR	pDONR223
NCBI SYMBOL	NKG7
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagctctgccgtccctgcccctgctggggggctccctgggctgatgtctgctgattgctttgagcaccgalltctgtttgagctgtgggtcccaccactcagctcactcgggctctgccaacagggcatggggacatcatatcaggctacatccacgtgacgcagacttc
agcattatggctgtctgtggccctgggtgctcctgagctctcctgctcctcctcactgttccccagggcaccgcccgtgtctcaaccaccgagcctttgctcagccatctccatgggtggtggccatggcggtgtacaccagcagcgggtgggaccagcctccacac
ccccagatccagaccttctcctggtcctctacctgggctgggtctcagctatcctctgctctgtacaggtgcctgagcctgggtgctcactgtggcggtccccgtcctggtatgaaacctgtga