

C9orf6-chromosome 9 open reading frame 6 Gene

Reference: PTXBC015795

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C9orf6; CG-8; Simiate; protein Simiate; protein FAM206A; family with sequence similarity 206 member A		

For research use only.

GENE DESCRIPTION	chromosome 9 open reading frame 6
GENE ID	54942
ORIGIN SPECIES	Human
ACCESSIONS	BC015795
VECTOR	pDONR223
NCBI SYMBOL	C9orf6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctaccgagccgaagccgagccggtggtgccttcgctcgtgacgatactcactcgtggtacaaacggatgtcaaaggaaatttggaggaccactgtatactacagcacttaaccgaatatgtgtcatcacattggcagaatctcatccagttcttcaaagtggaaaa
caattaaaagcatttctatcagatcagtaaccaactgtagcagacttcagaacaaggctctgggaaatttaagcgggggcacagtttctaacagagcttgacctctctgtaagatttactgtcagatggggaagaataatactgtgtctagtgtgttagagagcgtttgatggaagtgaatga
aaacatctccataagccatctattcttcaagaaaagccatctactgaaggctacattgcagttgtgttaccctaaatttgaagaaagtaaaagcataacagaagggttactgacacaaaaacaatatgaagaagtcattggtgaaacgcattatgccacaacagctacgtcatga