

RWDD1-RWD domain containing 1 Gene

Reference: PTXBC015802

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CGI-24; PTD013; RWD domain-containing protein 1; DRG family-regulatory protein 2; RWD domain containing 1		

For research use only.

GENE DESCRIPTION	RWD domain containing 1
GENE ID	51389
ORIGIN SPECIES	Human
ACCESSIONS	BC015802
VECTOR	pDONR223
NCBI SYMBOL	RWDD1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacagattacggcgaggagcagcgcaacgagctggagccctggagtcacatctacccctgactcctcagattatcagaaaatccaccagcttcacattactgtgacgtctgaggctggagaaaatgatgaaactgtccagactaccccaagttacatacagtgaaaaatacca
gatgaagctccccctttatgaaatattctccaggaaaatctagaagataatgatgtctcagacattttaaaattactagcattacaggctgaagaaaatcttggtatggtgatgatttttactctagtacagctgtgcaagaaaaatgaaatgaaatagtagatcagataaaactagaagagaaga
agaaaaagaacaaaaaagaagaagcagaagaagctgaaaagcaattattccatgggtactccagttacaattgagaatttcttaaaatggaaagccaagtttgatgcagaactcttggaaattaaaaagaaaaaggaatgaagagaagaacagcaggaaaaataaattagtgggaaa
caactatttgaacagatcataatcttgacacatctgatatccagttcttgaggatgctggaacaacgtggaggtagatgagctttgttccaagaaatggatgacttgagctggaggatgatgaagatgatccagactataatcctgctgacccagagagtgactcagctgactaa