



Reference: PTXBC015822

|                |                                                                                                                                                                                        |                     |                |
|----------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                                    | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                                   | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | GIG12; HEL110; HLF2; lactotransferrin; epididymis luminal protein 110; growth-inhibiting protein 12; kaliocin-1; lactoferricin; lactoferroxin; neutrophil lactoferrin; talalactoferrin |                     |                |

**For research use only.**

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | lactotransferrin                                                                                  |
| <b>GENE ID</b>          | 4057                                                                                              |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC015822                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | LTF                                                                                               |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

[illegible]