



Reference: PTXBC015823

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GIG12; HEL110; HLF2; lactotransferrin; epididymis luminal protein 110; growth-inhibiting protein 12; kaliocin-1; lactoferricin; lactoferroxin; neutrophil lactoferrin; talalactoferrin		

For research use only.

GENE DESCRIPTION	lactotransferrin
GENE ID	4057
ORIGIN SPECIES	Human
ACCESSIONS	BC015823
VECTOR	pDONR223
NCBI SYMBOL	LTF
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1