

NT5C3-5'-nucleotidase, cytosolic III Gene

Reference: PTXBC015856

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NT5C3; P5'N-1; P5N-1; PN-I; POMP; PSN1; UMPH; UMPH1; cN-III; hUMP1; p36; cytosolic 5'-nucleotidase 3A; 5'-nucleotidase, cytosolic III; 7-methylguanosine nucleotidase; 7-methylguanosine phosphate-specific 5'-nucleotidase; cytosolic 5'-nucleotidase 3; lupin		

For research use only.

GENE DESCRIPTION	5'-nucleotidase, cytosolic III
GENE ID	51251
ORIGIN SPECIES	Human
ACCESSIONS	BC015856
VECTOR	pDONR223
NCBI SYMBOL	NT5C3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatccagaattccagaaaagttcagttcgaatcaagaaccctacaagagtagaagaattatctgtgtcttatacaaggagagctgccaaactcagataataacggacttggatgacactcagtagatttcataaaaggaaaagatgcccaacatgtcataatcattgacaac
tftaagctggttacagatgaatgtagaaaaaagtftatgcaactaaaggaaaaatactacgctattgaagttgacctgttctactgtagaagagaaagtacccttatafgtggaatggtatacactaaatcacatggttctgttcagcgaagcttaccataaagctaaactaaagaaattgtggca
gaatcgcgttatgctcaagaagatgagaattctttgataagctccaacaacatagcatalcccggttccatatttccgctggaatcgccgatgtagaggaagttatcgtcaagctggtttatcccaatgtcaaaagtgtgtccaatttatgatttggatgaaactgggtgct
caaaggattaaaggagaactaattcatgtatttaacaacatgatgtgctctgaggaaatacagaatattcaatcaactaaaagacaatagtaacataattctctggagactccaaggagactaaagaatgacagatggagtgccaatgttgagcacattctgaaaattgatatcaaat
gatagagtgatgagcttttagaaaagatcatgactcttatgatattgttttagtacaagatgaatcattagaagtagccaactctattttacagaagattctataa