

CHRAC1-chromatin accessibility complex 1 Gene

Reference: PTXBC015891

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CHARC1; CHARC15; CHRAC-1; CHRAC-15; CHRAC15; YCL1; chromatin accessibility complex protein 1; DNA polymerase epsilon subunit p15; chromatin accessibility complex 15 kDa protein; histone-fold protein CHRAC15; chromatin accessibility complex 1		

For research use only.

GENE DESCRIPTION	chromatin accessibility complex 1
GENE ID	54108
ORIGIN SPECIES	Human
ACCESSIONS	BC015891
VECTOR	pDONR223
NCBI SYMBOL	CHRAC1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggacgtggtcgtgggttaagacaagggcggggagcagcgggtcatctcgtgcctctatcccgatccgggtcatcatgaagagctccccgaggtgtccagcatcaaccaggagcgttggtgctcacggccaaggccacggagctctttgttcaatgcctagccacctattc
ctacagacacggcagtggaagaaaaagaaagtactgacttacagtgtatttagcaaacactgcacagcaatcagaacttttcagtttcttgcagatataftaccaagaagatttagctagtaaatacctgaaatgcttaagaggaaaagggaagaagatgaggagaatgacaatga
taatgaaagtaccatgatgaagctgactcctaa