

ZNF684-zinc finger protein 684 Gene

Reference: PTXBC016016

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	zinc finger protein 684		
For research use only.			

GENE DESCRIPTION	zinc finger protein 684
GENE ID	127396
ORIGIN SPECIES	Human
ACCESSIONS	BC016016
VECTOR	pDONR223
NCBI SYMBOL	ZNF684
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaaaaatcacaaggaaaaaacctttgaatgcaatgactgtggaagcctatagcaggagggcacacttgcactcatcagaaaatcacaatggagagagacccttgtgtgcaatgattgtgggaaggcgttatgcataaagcccaactcgtggccaccagagacttcac
actggagagaagccttatgagtcagtcactgtgggaaacattcacttggaaactcctcatttaacacacgtgaaatctacacacttgagaagtcattgaaatgaaggaaatgtgggaaaccttcaggatagtcaccccttataaacattccagattcacaaggagagaacccctacc
agtgtatcatatgtggcaagcctttggcaacacacccgtgctgttacacaccaaagaattcacagagagaacccctacagttgtattgaaatgtggcaagccttcacagaagtcacatcctcagacatcagataactcacagagagaagccctatgaatgtaacagatgtggg
aaagcattttccagaagtcacaaatctattgtacatcagaaaattcacataaa