



Reference: PTXBC016041

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BP-4; HT29-IGFBP; IBP4; IGFBP-4; insulin-like growth factor-binding protein 4; IBP-4; IGF-binding protein 4; insulin like growth factor binding protein 4		

For research use only.

GENE DESCRIPTION	insulin-like growth factor binding protein 4
GENE ID	3487
ORIGIN SPECIES	Human
ACCESSIONS	BC016041
VECTOR	pDONR223
NCBI SYMBOL	IGFBP4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1