

CLDN3-claudin 3 Gene

Reference: PTXBC016056

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C7orf1; CPE-R2; CPETR2; HRVP1; RVP1; claudin-3; CPE-R 2; CPE-receptor 2; Clostridium perfringens enterotoxin receptor 2; ventral prostate.1 protein homolog; ventral prostate.1-like protein; claudin 3		

For research use only.

GENE DESCRIPTION	claudin 3
GENE ID	1365
ORIGIN SPECIES	Human
ACCESSIONS	BC016056
VECTOR	pDONR223
NCBI SYMBOL	CLDN3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtccatgggctggagatcagggcaccgcgtgcccgtgctgggctggctgggcaccatcgtgtgctgcgcgttgcctcgtggcgcgtgtgcgccttcacgcagcaacatcatcagctgcagaaacatctgggaggcctgtggatgaactgcgtggtgcagagcaccggcc
agatgcagtgcaaggtgtacgactcgtgctggcactgccacaggacgttcaggcggcccgccctcctcgtgtggccatcctgctggcgccttcgggctgtagtggcgctggtgggcgccagtgccacaaactgcgtgcaggacgacacggccaagccaagatcaccatcg
tggcaggcgtgctgttcttctcgcgcctgctcaccctcgtccgggtgctcgtggcggcaacacattatccgggacttctacaacccctggtgcccaggcgcagaaagcgcgagatgggcgcgggctgtacgtggctgggcggccggcgctgcagctgctggggggcg
cgctgctctgctgctgtgtccccacgcgagaagaagtacacggccaaggtcgtctactcgcgcgcgctccaccggccgggagccagcctgggcacaggtacgacgcaaggactacgtctaa