

## FAM120C-family with sequence similarity 120C Gene

Reference: PTXBC016138

|                |                                                                                                                                      |                     |                |
|----------------|--------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                  | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                 | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | CXorf17; ORF34; constitutive coactivator of PPAR-gamma-like protein 2; tumor antigen BJ-HCC-21; family with sequence similarity 120C |                     |                |

**For research use only.**

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | family with sequence similarity 120C                                                              |
| <b>GENE ID</b>          | 54954                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC016138                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | FAM120C                                                                                           |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atggacgccatgctgggctactgtcagcgtgtgccaggctgtgcctatctggcgccgacggcctggagctcgtgtcatgttccggggggcctgggcaaggaccggctggccgagtggtgggcccgtgtgccaggccgagcggcagacagcgaactgatcgtgggacacg  
tgggcaacaaggcaccctccaccgcgggctgttctgccaccggcctgagccactgcgtgaggttagcactatccgttccgggtcaagcaagctgccatgtgtgagctgccctatggagagggccacaggaaaagaactgagggaagcctccaaccaactcag  
tgaggaaactgagggccctcagtcctaactacagcctgcaaggaaactga