

CBY1-chibby homolog 1 (Drosophila) Gene

Reference: PTXBC016139

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C22orf2; CBY; HS508I15A; PGEA1; PIGEA-14; PIGEA14; arb1; protein chibby homolog 1; ARPP-binding protein; PKD2 interactor, Golgi and endoplasmic reticulum-associated 1; chibby CTNNB1-mediated transcription inhibitor; chibby homolog 1; coiled-coil protein P		

For research use only.

GENE DESCRIPTION	chibby homolog 1 (Drosophila)
GENE ID	25776
ORIGIN SPECIES	Human
ACCESSIONS	BC016139
VECTOR	pDONR223
NCBI SYMBOL	CBY1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccctttcttgggaatacgttcagtcggaagaagacacctcctcggagtcggcatctctccaactgcattcttggatcgaacccgggaggtggagctgggttggatacggatcccgactatgaacctggcaggcgcaaacctgaagttgaaaatggccagtggatagcagacaggggttagtggcggtgtggaccggagggttcagcgccttcgagcggaaccagcagttggaggagaagaacaatctcttgcggctgaagtggacatcttattagacatgctttcagagtcactgctgaatccacttaatggagaaggaaactggatgaactgaggatcagccggaagagaaaatga