

DLX4-distal-less homeobox 4 Gene

Reference: PTXBC016145

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BP1; DLX7; DLX8; DLX9; OFC15; homeobox protein DLX-4; beta protein 1; distal-less homeo box 7; distal-less homeo box 9; homeobox protein DLX-7; homeobox protein DLX-8; distal-less homeobox 4		

For research use only.

GENE DESCRIPTION	distal-less homeobox 4
GENE ID	1748
ORIGIN SPECIES	Human
ACCESSIONS	BC016145
VECTOR	pDONR223
NCBI SYMBOL	DLX4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaccttttgcctgccccctccccggcgggacgctcacaagctgtcttccagacctgccccctgtccctgtagcggtgcttaccctgcttgccttaccctacaccgagcagcgaaccccgagactctacctgtcttccagcaaccccgcgctctctcagccccctgtcggaacctgcagagcaccctcaggaactcagggcagactcggagaagccggtgtccccgggaacctccgagcgcgccctcaggcccccgccaaaagctccgcaagccgaggaccatctactccagcctgcagctgcagcacctaaaccagctttccagcacacgcagtaacctggcgctgcccagagggcccaagctggcagcgagctcggcctcaccagaccaggtaaagatctgtttcagaacaacgctccaagtataagaagctcctgaagcagaattctggggggcaggaaggggacttccctgggaggaccttctgtgtctccctgctccccaccctccctccctctgggatctaccaaggcaggaccctgcccaccagttggctatggcaacagctttggagcctggatcagcatcctcctcagatgtctggcttcgctcagatgatgtga