

C6orf136-chromosome 6 open reading frame 136 Gene

Reference: PTXBC016167

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C6orf136; chromosome 6 open reading frame 136		
For research use only.			

GENE DESCRIPTION	chromosome 6 open reading frame 136
GENE ID	221545
ORIGIN SPECIES	Human
ACCESSIONS	BC016167
VECTOR	pDONR223
NCBI SYMBOL	C6orf136
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaggaacatctgtctcatgtatgagagactgagacaagagcttcccaagctcttccttcagtcacgactacagctgtattccttggatgtggaaatcatcaatgagatcctcaacatactaccaaggccggacatggatcattcttctactgacctctgccgttctggcctgga
attattttgacacacctcgtttggaggttttacagctgacccgccaccctgagaactggacctgcaagcccggtggcgcttggggctgcccgtccacttgcctttttgcgggttctacaagcgtgacaaagacgagcattaccggacctatgatgcctactccacttctacctgaattccag
tggcctcatttgcgcatcgtctagataaactgatgccttcacactcacctccaacgcctgtgaagaagctgctagtggggagccctggtggccctggggctgtcagagccagaacctgacttaacctgtgttccaagccctga