

RPS2-ribosomal protein S2 Gene

Reference: PTXBC016178

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	LLREP3; 40S ribosomal protein S2; 40S ribosomal protein S4; OK/KNS-cl.6; protein LLRep3; ribosomal protein S2		

For research use only.

GENE DESCRIPTION	ribosomal protein S2
GENE ID	6187
ORIGIN SPECIES	Human
ACCESSIONS	BC016178
VECTOR	pDONR223
NCBI SYMBOL	RPS2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggatgacgcccgtgacgccccggggggccggggccctgtgtgcccctgggatgggaaccggtgtgctccgagggtttcggcagtgccatccggggccgggtcgcggccgtggacggggccggggccgaggccgcggagctcgcggaggcaaggccgagg
ataaggagtggaatgcccgtcaccaagttggccgcttggtcaaggacatgaagatcaagtcctggaggagatctatctctccctgccatlaaggatcagagatcattgatttctctggggccctctcctcaaggatgaggttttgaagattatgccagtcagaagcagacccgtgc
cggccagcgcaccaggttcaaggcattgtgtctatcggggactacaatggccacgtcgggtctgggtgttaagtgtcctcaaggaggtggccaccgccatccgtggggccatcatcctggccaaagctctccatcgtccctgtgcgagaggctactgggggaacaagatcggcaagcccc
acactgtcccttcaagggtgacaggccgtgctgtgtgtgacgctcatcctgcaccagggggactggcatgtctccgcacctgtgcttaagaagctgctcatgatggctggtatcgaactgtacacctcagccggggctgcactgccacctgggcaacttcgcaa
ggccacctttgatgccatttcaagacctacagctacgtgacccccgacctctggaaggagactgtattcaccaggtctccctatcaggagttcactgaccacctcgtcaagaccacaccagagctcctgtgcagcggactcaggtccagctgtggtacaacatag