

CD48-CD48 molecule Gene

Reference: PTXBC016182

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD48 molecule; CD48 antigen (B-cell membrane protein); CD48 antigen; BCM1; BLAST; BLAST1; MEM-102; SLAMF2; hCD48; mCD48; B-lymphocyte activation marker BLAST-1; BCM1 surface antigen; SLAM family member 2; TCT.1; leukocyte antigen MEM-102; signaling lympho		

For research use only.

GENE DESCRIPTION	CD48 molecule
GENE ID	962
ORIGIN SPECIES	Human
ACCESSIONS	BC016182
VECTOR	pDONR223
NCBI SYMBOL	CD48
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtgtccagaggttgggattcgtgtctgtgctctggaattgctactgtgctctgtcactcctggtgaccagcattcaaggtcacttggatcatatgacctgtgtctccggcagcaacgtgactctgaacatctctgagagcctgctgagaactacaacaactaacctggttttactttcg
accagaagattgtagaatgggattccagaaaatctaagtactttgaatccaaattaaaggcaggggtcagacttgatcctcagagtggtgcgactgtacatcttaaggtccagaaaggagacaacagcacctacatcatgaggggtgttgaataagactgggaatgagcaagaatggaagatc
aagctgcaagtgctgacctgtacccaagcctgtcatcaaaattgagaagatagaagacatggatgacaactgttattgaaactgtcatgtgtgatacctggcgagtgctgtaaaactacacgtgtatggggacaaaaggcccttccaaaggagctccagaacagtgctgtgaaaccacc
cttatgccacataattactccaggtgttatacttccaagtcagcaattctgtgagcagcaagaatggcaggtgtgctcagtcacacctgtaccctggcccggtccttggagtagaatggattgcaagttggctagtgtgtcacggtgccaccattctgacctgttacttactga