

CENPQ-centromere protein Q Gene

Reference: PTXBC016279

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C6orf139; CENP-Q; centromere protein Q		
For research use only.			

GENE DESCRIPTION	centromere protein Q
GENE ID	55166
ORIGIN SPECIES	Human
ACCESSIONS	BC016279
VECTOR	pDONR223
NCBI SYMBOL	CENPQ
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaagatttaactaatgtatcaagtctactgaatatggaaggcgcagacaaagctaatgaagaaggtctggcattactacaggaagaataagataaaatggtagagaccacagagttaatgactgggaatattcagagcctaaagaacaaattcagattctggcaagtgggtgga
agaagaagaggagagatgaacagatgcacaaataatagtagtgagtagtctctctccggaacttctcagaaaactctcaaaagcaccacactcagaaagaatttggcgctaatccaaaccagaatgctcttctaaaggactggatattctcataattcatcacagatgaagag
catgtcaacttcattgaagaagcctataagaactggatgcatttaa