

FKBP3-FK506 binding protein 3, 25kDa Gene

Reference: PTXBC016288

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PPIase FKBP3; peptidyl-prolyl cis-trans isomerase FKBP3; FKBP-25; FKBP-3; FKBP25; PPIase; 25 kDa FK506-binding protein; 25 kDa FKBP; FK506 binding protein 3, 25kDa; FK506-binding protein 25, T-cell; immunophilin FKBP25; rapamycin binding protein; rapamycin		

For research use only.

GENE DESCRIPTION	FK506 binding protein 3, 25kDa
GENE ID	2287
ORIGIN SPECIES	Human
ACCESSIONS	BC016288
VECTOR	pDONR223
NCBI SYMBOL	FKBP3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggcgccgtccacagcggcggtgacgctggagcagctgcgcagctgagcagctgcccagaaggacattcaagttctgcaggaacacgggtcagattcgttcttcgagaacataattattaggaacattaaaaatgtggccaagacagctaacaaggaccacttggtta
cagcctataccatcttttgaactaagcgttttaagggtactgaaagtataagtaagtgtctgagcaagtaaaaaatgtgaagcttaataagataaacccaaagaacaaagctggaagagaccctggatgagggtccacaaaataactaatctgttctgaaaaaggagataaaac
caacttcccaaaaaggagatgtgttcactgctgtatatacaggaactacaagatgggactgttttgataactaatattcaacaagtgcaagaagaagaaaaatgccaaagcctttaagttttaaggtcggagtaggcaagttatcagaggatgggatgaagctcttctgactatgagtaa
aggagaaaaggctcactggagattgaaccagaatgggttacggaagaaaggacagcctgatgccaaattccacaaatgcaaaactacttttgaagtgggaattagtgatattgattga